

Accurate, Handsfree Monitoring for Cardiac Stress and Exercise Testing



+ Reliable, automated blood pressure, + Dual blood pressure measurement so you can focus on your patient

modes, Exercise & Non-Exercise

SunTech Medical®'s Tango® M2 automated NIBP monitor uses DKA technology specifically engineered for the cardiac stress testing environment. With advanced proprietary algorithms SunTech Medical stands alone in the market of motion tolerant BP measurement. The Tango M2 allows for the focus to be on the patient by automatically taking accurate and reliable BP measurements rather than spending valuable time trying to take measurements manually. Tango M2 can add automated BP and SpO₂ to your stress ECG system, creating a complete cardiac diagnostic testing center.

+ Features & Benefits + Description

Seamless Integration: Non-Exercise BP Mode: Stat Mode: Color LCD: Verify Measurements: Data Retrieval: Field Updgrades: E-Library:

Automated communication with your stress system reduces the risk of transcription errors. Allows BP measurement during patient set-up and recovery without an ECG signal. Rapidly repeated automatic BP measurements for time sensitive and emergent situations. Improved usability with a new 7" color LCD display. See & hear the Korotkoff sounds using the onscreen display and the included headphones. Easier troubleshooting with 300 BP reading history and USB capabilities for measurement retrieval. USB port allows for field upgrades, making sure end-users always have the current software. Onboard E-Library, including step-by-step interface notes and tutorials.

+ Options

Single Patient Use (SPU) kits for increased infection control Pulse oximetry (SpO₂) Internal ECG



+ Specifications	+ Description
BP Measurement:	Auscultatory R-wave gating using K-sound analysis, for all static & active phases of stress testing.
	Oscillometric using pneumatic pressure for static measurements only.
Measurement Range:	Pressure (DKA Mode) Systolic: 40 - 270 mmHg, Diastolic: 20 - 160 mmHg (OSC Mode) Systolic:
	40 - 260 mmHg, Diastolic: 20 - 160 mmHg; Heart Rate: 40 - 200 bpm
Interfaces:	Integrates with all popular stress ECG systems using RS-232, BNC, ECG cable only for internal
	ECG & USB connections.
ECG Source:	Primary - From integrated stress ECG system or other external source
	Secondary - Internal ECG option using V2, V6, RL
Power:	Input - 100-240 VAC @ 1.5A, 50-60 Hz. Output: +9 VDC @ 5A IEC 320 type input connector
	Classification - Class II, continuous
BP Sampling Intervals:	From integrated stress ECG system or other external source, or 1 -20 minute intervals.
Dimensions:	24.0 cm x 17.4 cm x 11.5 cm (9.5" x 6.9" x 4.5")
Weight:	1.68 Kg (3.725 lb) (59.6 oz)
Warranty:	2 year standard warranty on monitor.
Accuracy:	Equivalent to a trained observer using a cuff/stethoscope auscultation method per ANSI/AAMI/ISO 81060-2
Standards:	IEC 60601-1:2005, IEC 60601-1-2:2007 (EMC), IEC 60601-2-30:2009, ISO 80601-2-61:2011, ISO
	10993-1:2009, ISO 10993-5,2009, ISO 10993-10:2010, FDA 21CFR801.5, MDD, WEEE

Autorizovaný obchodní a servisní partner:



COMPEK MEDICAL SERVICES, s.r.o. 17. listopadu 861, 506 01 Jičín nobil: + 420 605 281 433 el./fax: + 420 493 524 534 mail: info@compek.cz

COMPEK MEDICAL SERVICES, s.r.o. Strážna 11, 831 01 Bratislava mobil: + 421 908 758 793 fax: + 421 2 3301 6145 e-mail: info@compek.sk internet: www.compek.sk



507 Airport Boulevard, Suite 117 Morrisville, NC 27560-8200 USA Tel: + 1.919.654.2300 1.800.421.8626 Fax: + 1.919.654.2301

www.SunTechMed.com

SunTech Medical, Ltd. Oakfield Estate Eynsham, Oxfordshire England Tel: + 44 (0) 1865.884.234 Fax: + 44 (0) 1865.884.235

www.SunTechMed.com

SunTech Medical, Ltd.

www.CN.SunTechMed.com

2/F of Building A, Jinxiongda Technology Park Guanlan, Bao'an District Shenzhen, 518110, PRC Tel: + 86-755-2958 8810 + 86-755-29588986 (Sales) + 86-755-29588665 (Service) Fax: + 86-755-2958 8829

